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Precellys Plant DNA Kit (Safety-Line)

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CONTENT

INTRODUCTION	1
PRINCIPLE	1
KIT COMPONENTS	2
STORAGE AND STABILITY	2
BINDING CAPACITY	2
BEFORE STARTING	3
Initial preparation of Wash buffer	3
Required equipment and materials to be supplied by the user	3
DANGEROUS COMPONENTS	4
PRECELLYS PLANT DNA KIT PROTOCOL	6
CONCENTRATING THE DNA	8
QUANTIFICATION AND STORAGE OF THE DNA	8
TROUBLESHOOTING	9

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INTRODUCTION

The Precellys Plant DNA Kit provides a rapid and easy method for the isolation of up to 30 µg of genomic DNA from different tissues of a wide variety of plant species. The kit allows processing of single or multiple samples simultaneously with up to 100 mg of plant material per reaction.

The kit is specially developed for use with the homogenizer Precellys® 24 for a fast, effective and reproducible homogenization and cell disruption of a wide range of samples including hard to lyse tissues. Provided with the kit are 2 ml screw-cap tubes prefilled with a mixture of optimised amounts of ceramic beads of 1.4 mm and 2.8 mm in diameter giving best results for the lysis of plant tissues.

There is no need for extractions with organic solvents like phenol or chloroform or time-consuming steps such as precipitation with isopropyl alcohol or ethanol. Genomic DNA purified with the Precellys Plant DNA Kit is ready for a wide range of applications such as PCR, Southern Blotting and restriction digestion.

The Precellys Plant DNA Kit contains PerfectBind spin columns for the isolation of highly purified DNA. The columns can be closed tightly by lids to avoid cross-contamination more effectively.

PRINCIPLE

The Precellys Plant DNA Kit uses the reversible binding properties of the PerfectBind matrix, a silica-based material, combined with the speed of mini-column spin technology. A specifically formulated buffer system allows the isolation of up to 30 µg of genomic DNA with fragment sizes up to 60 kbp to bind to the matrix.

Samples are first homogenized with the homogenizer Precellys® 24, lysed under appropriate conditions and then applied to the PerfectBind® spin columns, where the DNA is effectively bound to the silica membrane. Cellular debris, proteins and other contaminants are washed away by specific buffers. The high quality DNA is finally eluted in Elution Buffer or sterile deionized water.

KIT COMPONENTS

Precellys Plant DNA Kit	5 Preparations	50 Preparations
Order No. Safety-Line	12-7531-00	12-7531-01
Components		
Precellys® Ceramic Kit 1.4/2.8 mm	5	50
PerfectBind® DNA columns	5	50
Microfilter	5	50
2 ml Collection Tubes	10	100
Lysis Buffer PL1-P	2 x 1.2 ml	25 ml
Lysis Buffer PL2-P	0.8 ml	7.5 ml
DNA Binding Buffer	1.4 ml	12.5 ml
DNA Wash Buffer	3.2 ml	32 ml
Elution Buffer (10 mM Tris-HCl, pH 9.0)	2 x 1 ml	18 ml
TE Buffer (10 mM Tris-HCL pH 8.0, 1 mM EDTA)	1.5 ml	10 ml
Instruction manual	1	1

STORAGE AND STABILITY

The components of the kit remain stable at an ambient temperature of 22 – 25 °C for at least 12 months after the date of purchase.

BINDING CAPACITY

Each PerfectBind® DNA column can bind approximately 30 µg of DNA. Using more than 100 mg of tissue is not recommended.

BEFORE STARTING

Please read the protocol carefully before the first use of the Precellys Plant DNA Kit and keep all required materials available before starting the preparation.

Initial preparation of Wash buffer

- ! DNA Wash Buffer is delivered as a concentrate and must be diluted with absolute ethanol before use:

5 Preparations	Mix 3.2 ml DNA Wash Buffer with 4.8 ml 100 % EtOH.
50 Preparations	Mix 32 ml DNA Wash Buffer with 48 ml 100 % EtOH.

Diluted DNA Wash Buffer should be stored at room temperature. If stored at lower temperatures, it should be heated to room temperature before use.

- ! All centrifugation steps have to be performed at 22 – 25 °C.

Required equipment and materials to be supplied by the user

- ! Precellys® 24 homogenizer
- ! Thermo Mixer or Digital Dry Bath for 2 ml tubes at 65°C
- ! Microcentrifuge for 2 ml tubes
- ! 100 % ethanol
- ! Sterile pipette tips and centrifuge tubes
- ! RNase A (20 mg/ml in 10 mM Tris-HCl, pH 8.0)

DANGEROUS COMPONENTS

Precellys Plant DNA Kit

Components	Signal word / symbols	Dangerous components	H and P statements
Precellys® Ceramic Kit 1.4/2.8 mm	-	-	-
PerfectBind DNA Columns	-	-	-
Microfilter	-	-	-
2 ml Collection Tubes	-	-	-
Lysis Buffer PL1-P	--- (not necessary)	sodium dodecyl sulphate, CAS: 151-21-3	--- (not necessary)
Lysis Buffer PL2-P	DANGER 	guanidinium thiocyanate 25-50%, CAS: 593-84-0	H302+H312+H332, H314, H412, P101, P102, P103, P260, P280, P303+P361+P353, P305+P351+P338, P405, P501
DNA Binding Buffer	DANGER 	propan-2-ol 50-100%, CAS: 67-63-0	H225, H319, H336, P101, P102, P103, P210, P261, P303+P361+P353, P305+P351+P338, P405, P501
DNA Wash Buffer	-	-	-
Elution Buffer (10 mM Tris-HCl, pH 9.0)	-	-	--
TE Buffer (10 mM Tris-HCL pH 8.0, 1 mM EDTA)	-	-	-

H and P statements	Descriptions
H225	Highly flammable liquid and vapour.
H302+H312+H332	Harmful if swallowed, in contact with skin or if inhaled.
H314	Causes severe skin burns and eye damage.
H319	Causes serious eye irritation.
H336	May cause drowsiness or dizziness.
H412	Harmful to aquatic life with long lasting effects.
P101	If medical advice is needed, have product container or label at hand.
P102	Keep out of reach of children.
P103	Read label before use.
P210	Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P260	Do not breathe dust/fume/gas/mist/vapours/spray.
P261	Avoid breathing dust/fume/gas/mist/vapours/spray.
P280	Wear protective gloves/protective clothing/eye protection/face protection.
P303+P361+P353	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P405	Store locked up.
P501	Dispose of contents/container to ...

PRECELlys PLANT DNA KIT PROTOCOL

1. Homogenization and lysis of plant tissues

Transfer up to 100 mg of plant material (cut hard to lyse tissue in pieces) into a tube of the Precellys® Ceramic Kit 1.4/2.8 mm and add 100 µl of TE buffer. Place the tube in the tube holder of the Precellys® 24 and homogenize 1 x 20 s at 5000 rpm.

For hard to lyse tissues an increase of homogenization energy might be necessary. Note that excessive homogenization can lead to shearing of the genomic DNA. Depending on the tissue the homogenization protocol has to be optimized by preliminary experiments.

Open the tube and add 400 µl Lysis Buffer PL1-P and 10 µl RNase A [20 mg/ml]. Vortex for 5 s and incubate for 15 min at 65 °C while shaking. If a shaking device is not available, mix the solutions 3 to 4 times by vortexing shortly.

Add 100 µl Lysis Buffer PL2-P and vortex for 5 s. Incubate the tube for 5 min on ice.

Centrifuge for 5 min at maximum speed.

Assemble a microfilter in a 2 ml collection tube. Transfer the clear sample (supernatant) carefully onto the microfilter and centrifuge for 1 min at 10000 x g (12 000 rpm in a standard microcentrifuge). Discard the microfilter and keep the filtrate.

2. Loading and Binding

Add 0.5 volume DNA Binding Buffer (in most cases 225 µl) to the filtrate and mix by pipetting up and down.

Assemble a PerfectBind DNA column in a 2 ml collection tube. Transfer up to 750 µl of the sample onto the column and centrifuge for 1 minute at 10000 x g. Discard the flow-through liquid and repeat loading the column with the remaining lysate if necessary. Discard the flow-through and keep the collection tube. Repeat centrifugation if the sample has not been centrifuged completely.

3. Washing

Add 650 µl DNA Wash Buffer (completed with ethanol!) to the column and centrifuge for 1 minute at 10000 x g. Discard the flow-through and keep the collection tube.

Add a second 650 µl of the DNA Wash Buffer on the column. Centrifuge the column with the collection tube for 1 minute at 10000 x g. Discard the flow-through and keep the collection tube.

4. Drying (Important step! Do not reduce centrifugation time!)

Place the column in the empty collection tube and dry the column completely by centrifugation for 2 min at 10000 x g.

5. Elution

Place the column into a sterile 1.5 ml microfuge tube and add 100 - 200 µl Elution buffer. Incubate tubes for 3 min at room temperature. To elute DNA from the column, centrifuge for 1 min at 6000 x g.

Optional: Repeat the elution with a second 100 µl Elution buffer.

A second elution increases the final yield, but can decrease the concentration of the isolated DNA. The final yield can be increased together with the concentration of the DNA, if the eluate is used for a second elution. But note that the yield is about 30 % decreased compared to a second elution with a fresh aliquot of Elution Buffer. The yield can be slightly improved by incubating the column at 70 °C instead of room temperature before centrifuging.

CONCENTRATING THE DNA

Genomic DNA purified with the Precellys Plant DNA Kit can be further concentrated if required. Add NaCl to an end concentration of 0.1 M and 2 sample volumes of absolute ethanol. Mix thoroughly by vortexing and incubate for 10 minutes at -20°C . Centrifuge for 15 minutes at $10000 \times g$ and discard the supernatant. Add 700 μl of 80 % ethanol and centrifuge for 2 minutes at $10000 \times g$. Discard the supernatant, air-dry the pellet for 2 minutes and dissolve the DNA in 20 μl of Elution buffer or sterile, deionised water.

QUANTIFICATION AND STORAGE OF THE DNA

To determine the concentration and the purity of a DNA containing solution, the absorbance of a 10- to 50-fold diluted aliquot is measured at 260 nm and 280 nm in a spectrophotometer. One A_{260} unity corresponds to 50 μg of DNA per ml. The concentration can be determined as follows:

DNA concentration ($\mu\text{g}/\text{ml}$) = absorbance₂₆₀ $\times$ 50 $\times$ dilution factor

The $A_{260/280}$ ratio of pure nucleic acids is 2.0. Generally genomic DNA isolated with the Precellys Plant DNA Kit shows ratios between 1.7 and 1.9 corresponding to a purity of 85 % to 95 %.

Genomic DNA isolated with the Precellys Plant DNA Kit can be stored in Elution Buffer, 10 mM Tris-HCl (pH 9.0) or sterile, deionised water for several years at -20°C . Avoid frequent freezing and thawing as this will shear the DNA and result in a reduction of molecule size.

TROUBLESHOOTING

Problem	Possible Cause	Suggestions
Clogged column	Incomplete lysis	Optimise homogenization parameter or extend incubation time of lysis with Lysis Buffer PL1-P.
	Sample too large	Divide the sample into several Precellys tubes to homogenize. Distribute lysate on more than one column.
	Sample too viscous	Dilute lysate with 10 mM Tris-HCl, pH 8.0. Distribute lysate on more than one column.
Low DNA yield	Incomplete lysis	See above.
	Poor elution	Pipette elution buffer onto the center of silica membrane. Repeat elution or increase elution volume. Incubate column at 70 °C for 5 minutes with Elution Buffer before centrifugation.
	Improper washing	Wash Buffer concentrate must be diluted with ethanol before first use.
	Clogged column	See above.
No DNA eluted	Poor binding of DNA to the PerfectBind column	Mix sample thoroughly with DNA Binding Buffer prior to loading onto PerfectBind column.
	No ethanol added to Wash Buffer concentrate	Dilute Wash Buffer with the indicated volume of absolute ethanol before first use.
Low $A_{260/280}$ ratio	Resin from the column present in the eluate	Avoid centrifugation at speeds higher than specified. The material can be removed from the eluate by centrifugation and will not interfere with enzymatic reactions.
	Low DNA yield	See above.
Problems with downstream applications	Salt carry-over	Ensure that the DNA Wash Buffer is at room temperature and does not contain any precipitates.
	Ethanol carry-over	Dry the columns after the washing step for at least 2 minutes by centrifugation at maximum speed.

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VWR International GmbH