Kurzanleitung – Short protocol

Homogenisator – Homogeniser Precellys® 24/Precellys® 24-Dual

INHALT / CONTENT

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CONTENT

WARRANTY	16
PACKAGING LIST	16
SAFETY INSTRUCTIONS	16
SYSTEM OVERVIEW	17
Technical Data	18
Combinable parameters/Refractory period	19
Instructions for the use of 7 ml tubes in the Precellys® 24-Dual	19
INSTALLATION	20
Installation	20
Turning on/Programming	20
Sample preparation	21
Loading the tube holder	22
Starting a run	22
Finishing a run	23
CLEANING AND DECONTAMINATION	23
TROUBLESHOOTING	24
APPLICATION TIPS	26
General user tips	26
Which lysis kit matches the starting material?	27
Further user tips	27
TECHNICAL SUPPORT AND ORDERING INFORMATION	28
Homogeniser Precellys® 24 and spare parts	28
Precellys® lysis kits	28
Precellys RNA isolation kits	29
Precellys DNA isolation kits	29

WARRANTY

PEQLAB guarantees that the product you have received has been thoroughly tested before shipment and meets its published specification. However, immediately upon arrival, please check carefully that the shipment is complete and has not been damaged in transit. For missing parts or to report any kind of damage, please contact PEQLAB (see 'TECHNICAL SUPPORT AND ORDERING INFORMATION'). Please retain all packaging materials until the delivery has been completely checked since this will speed up the return of goods if required and reduce environmental impact. Any form of returns, replacements or credit notes must be agreed in advance by PEQLAB.

For the Precellys® 24/Precellys® 24-Dual, PEQLAB guarantees a warranty period of 24 months if the product has been used solely according to the instruction manual and if not agreed differently. No liability is accepted for loss or damage arising from incorrect use. PEQLAB's liability is limited to the repair or replacement of the unit or refund of the purchase price at PEQLAB's discretion. PEQLAB is not liable for any consequential damages. Excluded from warranty are the following parts: fuses, toric joint, containment seal and indented plate. PEQLAB as well as the manufacturer BERTIN TECHNOLOGIES reserves the right to alter the technical specifications of the system without prior notice. This will enable us to implement developments as soon as they arise.

PACKAGING LIST

Unless otherwise agreed or written in the delivery note, the Precellys® 24/Precellys® 24-Dual contains the following components:

- One homogeniser Precellys® 24 or Precellys® 24-Dual including indented plate(s), toric joint and power cord (Precellys® 24-Dual is including one 7 ml and one 2 ml indented plate)
- One backup toric joint
- One backup fuse
- A User Manual in English and French language with all information about installation, operation and maintenance of the homogeniser Precellys® 24/Precellys® 24-Dual
- This short version of the User Manual in German and English language with essential information about installation and operation

SAFETY INSTRUCTIONS


- Please read the User Manual carefully before operating the Precellys® 24/Precellys® 24-Dual homogeniser.
- Disconnect the system from power supply as soon as it comes in contact with liquids.
- The system must always rest on its four feet.
- Do not exchange any components of the device and/or components of the software.
- Do not operate the unit when the casing is removed or considerably damaged.
- Do not open the lid before the sample holder comes to rest.
- Thermal expansion of liquid samples due to frictional heat might lead to tube leakage. Therefore please consider the information on maximal sample volumes given under 'Technical Data' and 'Instructions for usage of 7 ml tubes in the Precellys® 24-Dual'.
- Do not process samples containing explosive or highly flammable substances.

SYSTEM OVERVIEW

The Precellys® 24 homogeniser allows to simultaneously lyse up to 24 samples in 2 ml tubes. Using the Precellys® 24-Dual up to 6 samples in 7 ml tubes or up to 12 samples in 2 ml tubes can be grinded at once within a few seconds. The lysis principle is based on grinding the sample material with ceramic, steel or glass beads at very high speed in 2 ml or 7 ml screw cap tubes available as so called 'Precellys® lysis kits' (see 'TECHNICAL SUPPORT and ORDERING INFORMATION').

Therefore the tubes are put in the rotor-like tube holder of the Precellys® 24/Precellys® 24-Dual and are exposed to a powerful three dimensional precision movement (3-dimensional eight). Because of enormous power reserves, the choice of adequate lysis beads and the adaptation of the programmable parameter agitation speed, time and number of the agitation intervals, best homogenisation results can be obtained for a multitude of even challenging starting material and applications (see 'APPLICATIVE TIPS').

The patented vacuum closure system of the Precellys® 24/Precellys® 24-Dual allows for the fast and convenient loading of the tube holder and avoids the use of any screws or tools. Sensors in the locking latch and tube holder are permanently monitoring the locking state and will prevent or stop operation in case tubes or safety lid are not locked properly. Free accessibility and material consistence of the sample room which is sealed completely to the interior of the device is allowing for an effective decontamination of the system.

Front view, lid opened

Back view, lid closed

Technical Data

Electric properties		
Voltage range	240 V/50 Hz or 120 V/60 Hz (European or US standard)	
Power input	< 1 kVA	
Specification of fuses	5 x 20 – T 10 A – H 250 V	
Classification	Class I	
Physical properties		
Width	290 mm	
Height	370 mm at closed lid, 600 mm at open lid	
Depth	375 mm	
Weight	29.8 kg	
Generation of noise	< 70 dB	
Environmental conditions		
Temperature	15 - 30 °C	
Humidity	15 - 85 % HR	
Altitude (above sea level)	< 2000 m	
Programmable parameters		
Storable programs	5	
Agitation speed	Precellys® 24:	2500 rpm (Low Speed Program)* and 4000* - 6800 rpm, in steps of 100 rpm
	Precellys® 24-Dual:	2500 rpm (Low Speed Program)** and 4000** - 6500 rpm, in steps of 100 rpm
Duration of an agitation interval	5 - 90 s, in steps of 1 s (see 'combinable parameters...')	
Repetition of an agitation interval	1 - 3	
Pause between agitation intervals	5 - 120 s, in steps of 1 s (see 'combinable parameters...')	
Rest time between two consecutive runs (refractory period)	Precellys® 24:	2 or 5 min (see 'combinable parameters...'), disengageable***
	Precellys® 24-Dual:	none***
Sample capacity		
Sample number/run	Precellys® 24:	up to 24 x 2 ml tubes
	Precellys® 24-Dual:	up to 12 x 2 ml or up to 6 x 7 ml tubes
Maximal sample volume	0.35 ml in 0.5 ml tubes	
	1.4 ml in 2 ml tubes	
	5 ml in 7 ml tubes	

*from firmware version 4.1 (12/2011) on, until then 5000 – 6800 rpm, in steps of 100 rpm

**from firmware version 4.1 (12/2011) on, until then 5000 – 6500 rpm, in steps of 100 rpm

***from firmware version 4.1 (12/2011) on, until then refractory periods of 2 or 5 min (see 'Combinable parameters/Refractory period')
THE ACTUAL FIRMWARE VERSION OF THE SYSTEM IS DISPLAYED WHEN PRESSING 'CANCEL' FOR 10 SEC AFTER SWITCHING ON THE UNIT.

Combinable parameters/Refractory period

For safe operation the system is allowing the following combinations of program parameters.

Refractory periods are rest times in which a consecutive run cannot be started. They prevent the system from overheating. However, the risk of overheating only exists during extensive high throughput applications using Precellys®24. In addition system temperature is always monitored by a temperature sensor which will shut off the motor if the temperature rises above the limit. Therefore, from firmware version 4.1 (12/2011) on refractory periods do not exist any more in Precellys® 24-Dual and can be skipped by pressing 'cancel' in Precellys® 24.

The actual firmware version of the system is displayed when pressing 'cancel' for 10 sec after switching on the unit.

Speed [rpm]	Repetition of an agitation interval	Duration of an interval [s]	Pause between intervals [s]	Refractory period*
2500	1 – 3	5 – 90	5 – 120	2 min
4000 - 6100	1 – 3	5 – 90	5 – 120	2 min
6200	1 – 3	5 – 90	5 – 120	5 min
6300 - 6500	1 – 3	5 – 30	5 – 120	5 min
		31 – 60	15 – 120	
6600 – 6800*	1 – 3	5 – 30	30 – 120	5 min

* Precellys® 24 only

Instructions for the use of 7 ml tubes in the Precellys® 24-Dual

To prevent tube breakage caused by thermal expansion please comply to the following program parameters if using the 7 ml tubes in the Precellys® 24-Dual.

Speed [rpm]	Repetition of an agitation interval	Duration of an interval [s]
2500	1 – 3	5 – 60
4000 - 5500	1 – 3	5 – 40
5600 - 6000	1 – 3	5 – 30
6100 - 6500	1 – 2	5 – 20

INSTALLATION

Installation

Place the Precellys® 24/Precellys® 24-Dual onto a stable and level surface. Although the system has been manufactured in a most durable and solid way, resulting in a total weight of almost 30 kg, vibrations can occur during operation. Therefore please avoid to place the Precellys® 24/Precellys® 24-Dual close to devices/applications sensitive to any vibration (e.g. microscopy).

In order to provide optimal ventilation a minimum distance of 15 cm to the next device or to the wall should be kept.

Before connecting the power cord check that the voltage switch located at the back of the unit matches the voltage delivered by your electric power supply.

Place the delivered toric joint onto the tube holder. Be careful of the rubber ring not being porous or dirty.


Installing the toric joint into the notch of the tubes holder.


Turning on/Programming

Turn on the Precellys® 24/Precellys® 24-Dual by pressing the ON/OFF button located on the rear panel near the AC power input.

A welcome screen is displayed for a few seconds showing the device name Precellys® 24/Precellys® 24-Dual, the number of runs done by the device (R:) and the operating time in minutes (T:) corresponding to these runs.

Automatically the following basic menu is displayed showing program 1 and program 2 with the last defined parameters.


Move the cursor with the help of the arrow keys to the program of your choice and confirm with 'valid'.

The display now shows the selected program. Depending on the position of the double cursor you can start the program or change the parameter. By using 'cancel' the previous main menu appears.

After moving the cursor in the upper row and confirming with 'valid' the parameters of the program can be set. From the left to the right the parameters speed, number and duration of agitation intervals and the rest time between an agitation interval can be defined within the limits described in 'technical data'. For that purpose the flashing (= active) value has to be changed with the arrow keys and the set value confirmed with 'valid'. After that the next parameter starts to flash and can be set.

The programming is finished after having confirmed the value for the rest time of the intervals [PAUSE (s)]. The parameters of the changed program are shown in the display, the double cursor points at 'START', so that the program run can be directly started by confirming with 'valid'.


The flexible programming of the Precellys® 24 allows to adjust and optimise the lysis parameters for the effective lysis of various starting material. You will find some support in 'APPLICATIVE DETAILS'. We are also pleased to give you personal help.

Sample preparation


Put your sample/lysis buffer into a 2 ml or 7 ml screw cap tube of the appropriate Precellys lysis kit and close the tube properly. To prevent tube breakage caused by thermal expansion please pay attention that the volume of sample and buffer will not exceed 1.4 ml in 2 ml tubes and 5 ml in 7 ml tubes. Using the big ceramic or steel beads (CK28 or MK28) the volume should be at least 500 µl, because the power of the beads – especially at higher speed or longer agitation time – could lead to tube breakage (see also 'TROUBLESHOOTING' and 'APPLICATIVE TIPS').

Loading the tube holder

In order to load the tube holder with the 2 ml tubes the white indented plate has to be removed. After loading place the indented plate centred on the toric joint to hold down the tubes. The aglet provides for a correct positioning of the indented plate. After starting a run a vacuum system will press the indented plate onto the tubes which are fixed this way. To build up and keep the vacuum the toric joint has to be clean, intact and properly positioned.

Note for Precellys® 24-Dual: Use only one type of tubes (2 ml or 7 ml) at the same time. The indented plate 24 x 2 ml of the Precellys® 24 cannot be used in the Precellys® 24-Dual and vice versa.


Starting a run

Close the safety lid of Precellys® 24/Precellys® 24-Dual properly until the locking handle is locked.

Choose one of the programs from the main menu and confirm the double cursor pointing on 'START' with 'valid'.

As long as the vacuum between tube holder and indented plate is built up, the correlating message is shown in the display.


This phase takes a few seconds.

As soon as the vacuum is set the program starts according to the defined parameters. The display shows the remaining run time.


A run can be cancelled by pressing 'cancel' at any time.

Do not open the lid during the run. In case of opening the lid a security mechanism will stop the agitation movement within 4 seconds.

Finishing a run

At the end of a run, the movement of the tube holder is stopped and the unit displays the main menu. The vacuum between indented plate and tube holder is released about 5 seconds after the unit stops. After opening the security lid the indented plate can be taken off for the removal of the tubes.

Subsequent to a run on a Precellys® 24 the refractory period is displayed. Restarting a run is not possible during this time. The length of the refractory period is dependent on previous parameters of the program and prevents an overheating of the system. However it can be skipped by pressing 'cancel' and does not occur on a Precellys® 24-Dual. See 'Combinable parameters/Refractory periods' for details.

A run can be cancelled at any time by pushing 'cancel'.

Do not open the lid during the run. In case of opening the lid, a security mechanism will stop the agitation movement within 4 seconds.

CLEANING AND DECONTAMINATION

All parts of the casing can be cleaned with a damp cloth moistened with water, a mild cleaning agent or 70 % alcohol.

Besides 70 % of ethanol or standard polyalcohols you can also use 1.85 % bleaching chloride, 10 % formaldehyde or 1 % peracetic acid for disinfection. Completely remove the disinfectants after a residence time of 30 min with a damp cloth moistened with water.

Do not use strong solvents like acetone or concentrated bases (soda or caustic soda). Make sure that no liquid penetrates into the interior of the device.

TROUBLESHOOTING

In the following section you will find explanations and possible troubleshooting for the alarm messages shown on the display (this page) or other phenomena (next page).

Alarm message	Possible cause	Suggestions
Pressure alarm	Tubes are inapplicable or not positioned properly.	Use the tubes of the Precellys lysis kits offered by PEQLAB. Reposition the tubes on the holder.
	The indented plate is not centred or placed even.	Reposition the indented plate.
	The toric joint is not positioned or positioned improperly.	Reposition the toric joint.
	The indented plate is damaged.	Replace the indented plate.
	The toric joint is dirty or damaged.	Clean the toric joint or replace it.
	The vacuum is still not established properly.	Turn off the system and contact the TECHNICAL SUPPORT at PEQLAB.
Top alarm	The security cover is not locked properly.	Check that nothing prevents the cover from closing. Press cover and ensure that the handle is locked properly.
	The detection system is faulty.	Turn off the system and contact the TECHNICAL SUPPORT at PEQLAB.
Err speed sensor	Speed regulation system or detection system is faulty.	Turn off the system and contact the technical service at PEQLAB.
Over temp wait!	Engine temperature has reached the safety limit.	Do not switch off the system so that the motor ventilation can work. Ensure air openings and the cooling fan are not blocked up. If the alarm is still on after 30 min of cooling, contact the TECHNICAL SUPPORT at PEQLAB.
	The temperature detection system is faulty.	Turn off the system and contact the TECHNICAL SUPPORT at PEQLAB.
Err power supply	Power supply is not suitable and wrong setting of the power input respectively.	Make sure that the power input filter on the back of the device matches the voltage range delivered by your local mains supply.
	No synchronization between unit and power supply.	Turn off the system and contact TECHNICAL SUPPORT at PEQLAB.

Problem	Possible cause	Suggestions
Impaired running smoothness, increased vibration	The four feet do not rest on a level surface.	Adjust the feet with the help of a screw wrench. Please work with someone together to carefully tilt the Precellys® 24/Precellys® 24-Dual. Never set the Precellys® 24/Precellys® 24-Dual aside because the weight could damage the casing.
The fan does not work	Electric power supply is disconnected.	Check the unit is plugged in properly. Make sure that the power input filter on the back of the device matches the voltage range delivered by your local mains supply.
	Faulty fuse	Replace fuse.
	Faulty cooling system	Turn off the system and contact TECHNICAL SUPPORT at PEQLAB.
No display on the screen	Electric power supply is disconnected	Check the unit is plugged in properly. Make sure that the power input filter on the back of the device matches the voltage range delivered by your local mains supply.
	Faulty fuse	Replace fuse.
	Faulty display	Turn off the system and contact TECHNICAL SUPPORT at PEQLAB.
One or several tubes are not tight	The cap is not screwed properly.	Lock the tubes properly. Clean and/or disinfect the system as described in 'CLEANING AND DECONTAMINATION'.
	Sample material or beads sitting between cap and tube prevent tight closing.	Occurs eventually using very small beads (e.g. VK01). Pay attention that there is no sample material or beads located at appropriate sites of the cap or tube (optionally tubes can be briefly centrifuged).
	Tubes are disrupted by the kinetic energy of the beads.	Use exclusively high-quality screw cap tubes e.g. the Precellys lysis kits of PEQLAB. By using the larger beads (e. g. CK28 or MK28) the sample volume should be at minimum 500 µl (see also 'General user tips')!
	Tubes are disrupted by gas pressure.	The maximum sample volume of 1.4 ml in the 2 ml tubes and 5 ml in the 7 ml tubes should not be exceeded in order to guarantee sufficient compressible volume for a potential heat expansion in the tube.
Power failure, pulling the main plug during the run.	For safety purposes the vacuum between indented plate and sample holder is held.	Make sure the electrical power supply is working and turn the system back on. The vacuum is relieved then and the indented plate can be removed.

APPLICATION TIPS

For the homogenisation of a wide variety of starting materials PEQLAB offers 19 different optimised Precellys lysis kits, each consisting of high-quality 0.5 ml, 2 ml or 7 ml screw cap tubes filled with ceramic, steel or glass beads. They are suitable for use in the MINILYS, Precellys® 24/Precellys® 24-Dual and other 'bead beating' homogenisers. The kits are also available as part of the Precellys RNA and DNA isolation kits designed for homogenisation and downstream isolation of RNA and DNA.

19 different Precellys® lysis kits available

General user tips

Essential for the selection of the appropriate lysis kit is the consistence of the starting material to homogenise.

Mechanically highly resistant samples e.g. hard and elastic samples like bones or cartilage require high lysis energy which can be generated by the larger ceramic or steel beads. For the 2 ml tubes they are available in two types 'standard' and 'reinforced': 1) when using transparent 'standard tubes' due to the enormous kinetic energy of the beads, to prevent tube damage, do not use less than 500 µl buffer. 2) the recently introduced 'reinforced tubes' are suitable for use with smaller sample volumes even at very high energy conditions. The use of liquid nitrogen will make the tubes brittle and is not recommended.

If the lysis buffers contain potentially corroding material (including phenol, chaotropic salts, detergents etc.) use ceramic beads instead of steel beads.

For most animal and plant tissue the 1.4 mm ceramic beads (CK14) or the 1.4/2.8 mm ceramic mix (CKM) are more suited. Although the beads have a large grinding potential damage to the reaction tubes does not occur. Additionally, the ceramic material is chemically extremely inert and therefore suitable for most lysis buffers and downstream applications.

The lysis kit SK38 with its optimized mixture of glass and ceramic beads is available for the effective lysis of microorganisms in soil samples.

For the disruption of microorganisms glass or ceramic beads with a diameter of 0.5 mm or 0.1 mm are ideal. As a rule of thumb 'The smaller the material to lyse the smaller the beads to use'.

For very small amounts of sample (e.g. weight < 10 mg or volume < 100 µl) the 1.4 mm ceramic beads or the 0.5 mm glass beads are available in smaller 0.5 ml tubes (CK14S or VK05S).

The following table provides an overview for choosing the appropriate lysis kit.

Which lysis kit matches the starting material?

	human & animal tissue								plant tissue					microorganisms/cells					
	Skeletal muscle	Lung, heart	Liver, kidney	Brain	Connective tissue	Hair	Finger nails	Bones	Crop seeds	Small seeds	Stem	Roots	Leaves	Soil samples	Fungi	Cell culture	Yeast	Bacteria	Spores
Ceramic Kit 5.0/2.8 mm, 7 ml (CK50L)					x	x	x	x	x	x	x	x							
Ceramic Kit 2.8 mm, 7 ml (CK28L)*	x	x	x	x	x	x	x	x	x	x	x	x							
Keramik-Kit 1.4 mm, 7 ml (CK14L)	x	x	x	x								x		x					
Glass/Ceramic Kit, 7 ml (SKL)*			x	x								x	x*	x					
Ceramic Kit 0.1 mm, 7 ml (CK01L)														x	x	x	x	x	x
Glass Kit 0.1/0.5 mm, 7 ml (VKML)														x	x	x	x	x	x
*if needed, grinding efficiency can be improved by adding single ceramic beads 5.0 mm or 6.8 mm (CK50P or CK68P) manually.																			
Steel Kit 2.8 mm, 2 ml (MK28, MK28R*)						x	x	x	x										
Ceramic Kit 6.8 mm, 2 ml (CK68R*)					x		x	x	x										
Ceramic Kit 5.0/2.8 mm, 2 ml (CK50R*)					x	x	x	x	x	x	x	x							
Ceramic Kit 2.8 mm, 2 ml (CK28, CK28R*)	x	x			x	x	x	x	x	x	x	x							
Ceramic Kit 1.4/2.8 mm, 2 ml (CKM)	x	x	x	x	x	x				x	x	x	x		x				
Ceramic Kit 1.4 mm, 2 ml (CK14)	x	x	x	x									x		x				
Glass/Ceramic Kit, 2 ml (SK38)														x	x	x	x	x	x
Glass Kit 0.5 mm, 2 ml (VK05)															x	x	x	x	x
Glass Kit 0.1 mm, 2 ml (VK01)																x	x	x	x
*MK28R & CK28R tubes are mechanically reinforced, however, prone to stronger embrittlement in liquid nitrogen.																			
Ceramic Kit 1.4 mm, 0.5 ml (CK14S)	x	x	x	x	x					x	x	x	x		x				
Glass Kit 0.5 mm, 0.5 ml (VK05S)															x	x	x	x	x
0.5 ml tubes are ideally suited for small sample volumes < 100 µl or < 10 mg.																			

Further user tips

To preserve the integrity of the target molecule homogenisation should never be more powerful than required (see Fig. A). For example, excessive grinding can cause shearing of genomic DNA (see Fig. B). Frictional heat occurring during homogenising generally has no negative effect on nucleic acids. However, any heat should be avoided if native proteins will be analysed. For this purpose the optional use of the **Cryolys Cooling Module** for the Precellys® 24/Precellys® 24-Dual is recommended. If using foaming lysis buffers grinding energy must be increased, because the beads are slowed down by the foam. The reagent peqGOLD TriFast™ has a very high lysis potential for the isolation of RNA, DNA or proteins. While its ingredients (amongst others phenol and GTC) effectively support cell lysis, it is almost non-foaming. For starting material like soil samples, bacteria, yeast, fungi, animal and plant tissue optimised kit combinations are available for the homogenisation and downstream nucleic acid isolation using silica-based centrifugation columns (see 'TECHNICAL SUPPORT and ORDERING INFORMATION').



Fig. B) Mechanical stability of DNA. Genomic DNA from *Arabidopsis* leaves was isolated using the peqGOLD Plant DNA Mini Kit (homogenised in the Precellys®24). Afterwards aliquots of the same amount were exposed increasing grinding energies in CK14 tubes using the homogeniser Precellys® 24 and separated on an agarose gel ('M' = peqGOLD 1 kb DNA-ladder).

TECHNICAL SUPPORT AND ORDERING INFORMATION

For technical questions please contact us by phone +49 (0)9131 610 7020 or e-mail service@peqlab.de. For detailed information on our products please visit our website www.peqlab.com.

Homogeniser Precellys® 24 and spare parts

Product	Description	Order no.
Homogeniser Precellys® 24	powerful homogeniser for the simultaneous lysis of up to 24 samples in 2.0 ml screw cap tubes	91-PCS24
Homogeniser Precellys® 24-Dual	powerful homogeniser for the simultaneous lysis of up to 6 samples in 7 ml or up to 24 samples in 2 ml tubes	91-PCS24D
Cryolys	cooling module for the Precellys® 24/Precellys® 24-Dual	91-PCS-CRYO
Indented plate for Precellys® 24	for the reliable fixation of the tubes by vacuum	91-PCS24-AP
Indented plate 2 ml (Precellys® 24-Dual)	for the reliable fixation of the 2 ml tubes by vacuum	91-PCS24D-A2
Indented plate 7 ml (Precellys® 24-Dual)	for the reliable fixation of the 7 ml tubes by vacuum	91-PCS24D-A7
Toric joint	for the airtight sealing of the indented plate	91-PCS24-DR
Containment seal	for sealing up the sample room to interior of the device	91-PCS24-GM
Fuse	fuse of the power input	91-PCS24-SI
Anti rotation kit	vacuum tubing with metal protection coating	91-PCS24-AK

Precellys® lysis kits

Product	Description	Order no.
Precellys steel kit 2.8 mm	50 x 2.0 ml tubes with steel beads, Ø 2.8 mm	91-PCS-MK28
Precellys steel kit 2.8 mm, reinforced*	50 x 2.0 ml tubes with steel beads, Ø 2.8 mm Highly impact resistant tubes*	91-PCS-MK28R
Precellys Ceramic kit 6.8 mm, reinforced*	50 x 2.0 ml Tubes with ceramic beads, Ø 6.8 mm Highly impact resistant tubes*	91-PCS-CK68R
Precellys Ceramik kit 5.0/2.8 mm, reinforced*	50 x 2.0 ml Tubes with ceramic beads, Ø 5.0 and 6.8 mm, Highly impact resistant tubes*	91-PCS-CK50R
Precellys ceramic kit 2.8 mm	50 x 2.0 ml tubes with ceramic beads, Ø 2.8 mm	91-PCS-CK28
Precellys ceramic kit 2.8 mm, reinforced*	50 x 2.0 ml tubes with ceramic beads, Ø 2.8 mm Highly impact resistant tubes*	91-PCS-CK28R
Precellys ceramic kit 1.4/2.8 mm	50 x 2.0 ml tubes with ceramic beads, Ø 1.4 and 2.8 mm	91-PCS-CKM
Precellys ceramic kit 1.4 mm	50 x 2.0 ml tubes with ceramic beads, Ø 1.4 mm	91-PCS-CK14
Precellys ceramic kit 1.4 mm 'small'	50 x 0.5 ml tubes with ceramic beads, Ø 1.4 mm	91-PCS-CK14S
Precellys glass/ceramic kit SK38	50 x 2.0 ml tubes with glass & ceramic beads	91-CKS-SK38
Precellys glass kit 0.5 mm	50 x 2.0 ml tubes with glass beads, Ø 0.5 mm	91-PCS-VK05
Precellys glass kit 0.5 mm 'small'	50 x 0.5 ml tubes with glass beads, Ø 0.5 mm	91-PCS-VK05S
Precellys glass kit 0.1 mm	50 x 2.0 ml tubes with glass beads, Ø 0.1 mm	91-PCS-VK01

Precellys ceramic kit 5.0/2.8 mm, 7 ml	50 x 7.0 ml tubes with ceramic beads, ∅ 5.0 & 2.8 mm	91-PCS-CK50L
Precellys ceramic kit 2.8 mm, 7 ml	50 x 7.0 ml tubes with ceramic beads, ∅ 2.8 mm	91-PCS-CK28L
Precellys ceramic kit 1.4 mm, 7 ml	50 x 7.0 ml tubes with ceramic beads, ∅ 1.4 mm	91-PCS-CK14L
Precellys ceramic kit 0.1 mm, 7 ml	50 x 7.0 ml tubes with ceramic beads, ∅ 0.1 mm	91-PCS-CK01L
Precellys glass/ceramic kit 0.1/1.4 mm, 7 ml	50 x 7.0 ml tubes with glass & ceramic beads, ∅ 0.1/1.4 mm	91-PCS-SKL
Precellys glass kit 0.1/0.5 mm, 7 ml	50 x 7.0 ml tubes with glass beads, ∅ 0.1/0.5 mm	91-PCS-VKML
Precellys ceramic beads 5.0 mm	200 pieces/ca. 80 g, ∅ 5.0 mm	91-PCS-CK50P
Precellys ceramic beads 6.8 mm	200 pieces/ca. 170 g, ∅ 6.8 mm	91-PCS-CK68P

*not suitable for the use of liquid nitrogen (increased brittleness)

Precellys RNA isolation kits

Article	Description	Order-No. (50 Preps)*
Precellys Bacterial/Fungal RNA Kit	for the isolation of up to 100 µg total RNA from bacteria, yeast, fungi and spores. Including 50 x 2.0 ml tubes with glass beads (∅ 0.5 mm).	12-7611-01
Precellys Tissue RNA Kit	for the isolation of up to 100 µg total RNA from animal tissue. Including 50 x 2.0 ml tubes with ceramic beads (∅ 1.4 and 2.8 mm).	12-7621-01
Precellys Plant RNA Kit	for the isolation of up to 100 µg total RNA from plant tissue. Including 50 x 2.0 ml tubes with ceramic beads (∅ 1.4 and 2.8 mm).	12-7631-01

*additional package sizes for 5 or 200 preparations available

Precellys DNA isolation kits

Article	Description	Order-No. (50 Preps)*
Precellys Soil DNA Kit	for the isolation of up to 30 µg DNA from soil samples. Including 50 x 2.0 ml tubes with glass & ceramic beads (SK38).	12-7501-01
Precellys Bacterial/Fungal DNA Kit	for the isolation of up to 30 µg DNA from bacteria, yeast, fungi and spores. Including 50 x 2.0 ml tubes with glass beads (∅ 0.5 mm).	12-7511-01
Precellys Tissue DNA Kit	for the isolation of up to 50 µg DNA from animal tissue. Including 50 x 2.0 ml tubes with ceramic beads (∅ 1.4 and 2.8 mm).	12-7521-01
Precellys Plant DNA Kit	for the isolation of up to 30 µg DNA from plant tissue. Including 50 x 2.0 ml tubes with ceramic beads (∅ 1.4 and 2.8 mm).	12-7531-01

* additional package sizes for 5 or 200 preparations available

D PEQLAB Biotechnologie GmbH, 91052 Erlangen, Freecall (D): 0800 100 20 16, info@peqlab.de, www.peqlab.de
AT PEQLAB Biotechnologie GmbH, 6404 Polling, Tel: +43 (0) 5238 84 169, info@peqlab.at, www.peqlab.at
UK PEQLAB Ltd., Southampton SO31 7ZN, Freephone (UK): 0808 202 1302, info@peqlab.co.uk, www.peqlab.co.uk
USA PEQLAB LLC, Wilmington, DE 19810, Toll-Free (US): 877 737 5220, info@peqlab.us, www.peqlab.us

Creating the future together.