

## Datasheet

### AGR2 monoclonal antibody (M01), clone 1E5

**Catalog Number:** H00010551-M01

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody  
raised against a full-length recombinant AGR2.

**Clone Name:** 1E5

**Immunogen:** AGR2 (AAH15503.1, 1 a.a. ~ 175 a.a)  
full-length recombinant protein with GST tag. MW of the  
GST tag alone is 26 KDa.

**Sequence:**

MEKIPVSAFLLLVALSYTLARDTTVKPGAKKDTKDSRP  
KLPQTL SRGWGDQLIWTQTYEEALYKSKTSNKPLMIH  
HLDECPHSQALKKVF AENKEIQKLAEQFVLLNLVYETT  
DKHLSPDGQYVPRIMFVDP SLTVRADITGRYSNRLYAY  
EPADTALLLDNMKKALKLLKTEL

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA, S-ELISA, WB-Ce, WB-Re, WB-Ti,  
WB-Tr

(See our web site product page for detailed applications  
information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product  
page for detailed protocols

**Isotype:** IgG2b Kappa

**Storage Buffer:** In 1x PBS, pH 7.4

**Storage Instruction:** Store at -20°C or lower. Aliquot to  
avoid repeated freezing and thawing.

**Entrez GeneID:** 10551

**Gene Symbol:** AGR2

**Gene Alias:** AG2, GOB-4, HAG-2, XAG-2