

Datasheet

HTATIP2 polyclonal antibody (A01)

Catalog Number: H00010553-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length recombinant HTATIP2.

Immunogen: HTATIP2 (AAH02439, 1 a.a. ~ 242 a.a)
full-length recombinant protein with GST tag.

Sequence:

MAETEALSKLREDFRMQNKSVFILGASGETGRVLLKEI
LEQGLFSKVTLIGRRKLTfDEEAYKNVNQEVVDFEKLD
DYASAFQGHdVGFCClGTTRGKAGAEgFVRVDRDYV
LKSAELAKAGGCKHFNLLSSKGADKSSNfLYLQVKGE
VEAKVEELKFDRYSVFRPGVLLCDRQESRPGEWLVRK
FFGSLPDSWARGHSVPVVTVVRAMLNNVVRPRDKQM
ELLENKAIHDLGKAHGSLKP

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re, WB-Ti
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 10553

Gene Symbol: HTATIP2

Gene Alias: CC3, FLJ26963, SDR44U1, TIP30