

Datasheet

DPYSL4 purified MaxPab mouse polyclonal antibody (B01P)

Gene Symbol: DPYSL4

Gene Alias: CRMP3, DRP-4, ULIP4

Catalog Number: H00010570-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human DPYSL4 protein.

Immunogen: DPYSL4 (NP_006417.1, 1 a.a. ~ 572 a.a) full-length human protein.

Sequence:

MSFQGKKSIPIRITSDRLLIRGGRIVNDDQSFYADVHVE
DGLIKQIGENLIVPGGIKTIDAHGLMVLPGGVDVHTRLQ
MPVLGMTPADDFCQGTKAALAGGTTMILDHVFPDTGV
SLLAAEQWRERADSAACCDYSLHVDITRWHESIKEE
LEALVKEKGVNSFLVFMAYKDRCQCSDSQMYEIFSIR
DLGALAQVHAENGDIVVEEQKRLELGITGPEGHVLSH
PEEVEAEAVYRAVTIAKQANCPYVTKVMSKGAADAIA
QAKRRGVVVFGEPIASLGTGSHYWSKNWAKAAAF
VTSPPVNPDPTTADHLTCLSSGDLQVTGSAHCTFTTA
QKAVGKDNFALIPEGTNGIEERMSMVWEKCVASGKM
DENEFAVTSTNAAKIFNFYPRKGRVAVGSDADLVIWN
PKATKIISAKTHNLNVEYNIFEGVECRGAPAVVISQGRV
ALEDGKMFVTPGAGRFVPRKTFPDFVYKRIKARNRLA
EIHGVPRGLYDGPVHEVMVPAKPGSGAPARASCPGKI
SVPPVRNLHQSGFSLSGSQADDHIARRTAQKIMAPPG
GRSNITSLS

Host: Mouse

Reactivity: Human, Rat

Applications: Det Ab, WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 10570