

## Datasheet

### ADAP1 mouse monoclonal antibody (hybridoma)

**Catalog Number:** H00011033-M

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody  
raised against a full-length recombinant ADAP1.

**Immunogen:** ADAP1 (NP\_006860.1, 1 a.a. ~ 374 a.a)  
full-length recombinant protein with GST tag. MW of the  
GST tag alone is 26 KDa.

**Sequence:**

MAKERRRAVLELLQRPGNARCADCGAPDPDWASYTL  
GVFICLSCSGIHRNIPQVSKVKSRLDAWEEAQVEFMA  
SHGNDAAARARFESKVPsfYYRPTPSDCQLLREQWIRA  
KYERQEFYIPEKQEPYSAGYREGFLWKRGRDNGQFL  
SRKFVLTEREGALKYFNRNDAKEPKAVMKIEHLNATFQ  
PAKIGHPHGLQVTYLKDNSTRNIFIYHEDGKEIVDFWN  
ALRAARFHYLQVAFPGASDADLVPKLSRNYLKEGYME  
KTGPKQTEGFRKRWFTMDDRRLMYFKDPLDAFARGE  
VFIGSKESGYTVLHGFPSTQGHHPHGITIVTPDRKF  
LFACETESDQREWVAAFQKAVDRPMLPQEYAVEAHF  
KHKP

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA, WB-Re, WB-Tr  
(See our web site product page for detailed applications  
information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product  
page for detailed protocols

**Entrez GeneID:** 11033

**Gene Symbol:** ADAP1

**Gene Alias:** CENTA1, GCS1L, p42IP4