

Datasheet

FDCSP polyclonal antibody

Catalog Number: PAB20711

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against recombinant FDCSP.

Immunogen: Recombinant protein corresponding to amino acids of human FDCSP.

Sequence:

KRSISDSDELASGFFVFPYPYPRPLPPIPFPRFPWFR
RNFPIPIESAPTTPL

Host: Rabbit

Reactivity: Human

Applications: IHC-P

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgG

Recommend Usage: Immunohistochemistry

(1:50-1:200)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 260436

Gene Symbol: FDCSP

Gene Alias: UNQ733/PRO1419, C4orf7, FDC-SP

Gene Summary: FDCSP is expressed by follicular dendritic cells (FDCs) and activated leukocytes during ongoing immune responses.[supplied by OMIM]