



ZPACK13260004 CUBITAINER 10 L UN



SPECIFICATIONS

UN-certification : 6HG2/Y1.4/Z1.4/30/.....

Material : PE

Filling Opening : 34,50 mm

Height : 232,00 mm

Length : 236,00 mm

Wide : 236,00 mm

Remark : Packaging materials and specifications shown may be modified without preceding notice
Version dd. 22 mai 2012

tel/tél +32(0)16 385 011
fax +32(0)16 385 385
info@be.vwr.com
www.vwr.com

Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001 Leuven
RPR/RPM Leuven
BTW/TVA BE 0403.593.343