



Product Specification Sheet

VWR FILTER TIPS BULK

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Type	Capacity (µl)	Length(mm)	Inner Diameter (mm)	Packed	Pk
732-0565	UltraFine™, extended	0,1 - 10 µl	38.1	0.25	Bulk	1
732-0561	UltraFine™, graduated	0,1 - 10 µl	31.4	0.51	Bulk	1
732-0517	Bevel Point™, graduated	0,5 - 10 µl	44.8	0.51	Bulk	1
732-0526	Bevel Point™	1 - 40 µl	49.6	0.51	Bulk	1
732-0567	Bevel Point™	1 - 50 µl	49.6	0.55	Bulk	1
732-0521	Bevel Point™	1 - 100 µl	49.6	0.55	Bulk	1
732-0571	Bevel Point™	1 - 160 µl	49.6	0.51	Bulk	1
732-0558	Bevel Point™, wide orifice	1 - 200 µl	49.6	1.45	Bulk	1
732-0539	UltraFine™	1 - 200 µl	52	0.51	Bulk	1
732-0535	Bevel Point™	1 - 300 µl	52	0.51	Bulk	1
732-0532	UltraFine™	100 - 1000 µl	88.9	0.76	Bulk	1

Material:

- Polypropylene made of 100% virgin resin
- The resins used to produce these products have been tested for heavy metals contamination using the U.S.P. (231) method and found to have results lower than 1ppm.



Product Specification Sheet

Other Characteristics:

- Made in the highest quality moulds, these tips can be used whenever the work demands a higher standard.
- A wider seal assures that tips will seal easily and consistently without leaking, on old as well as new pipettes
- Autoclavable

Clean Claims:

- Certified testing is performed lot by lot.
- Free of detectable Nuclease contaminants(RNase,DNase) as confirmed by a gel assay
- Free from Human DNA to a detection level of 1 pico-gram of DNA
- Certificates are on file at the production facility.

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **January 2014**