



Certificate of Conformity

VWR PP CARBOYS, NARROW MOUTH, ROUND, WITH SCREW CAPS, NON-STERILE

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Catalog No.	Capacity (L)	Neck I-Ø (mm)	Thread (mm)	General Wall Thickness (mm)	Product Size (±2mm)	Pk
216-1703	10	65	83	2.5	Ø250 X 385	1
216-1701	20	65	83	2.5	Ø287 X 535	1
216-1702	50	65	83	2.5	Ø378 X 680	1

Material

Carboy : Polypropylene, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed.

Closure : Polypropylene, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed With Food grade colour.

BSE/TSE Statement: These product is manufactured without the use of Raw materials of animal origin.

Latex Free Statement: No material containing latex or natural rubber are used in the manufacturing, handling and packaging process for this product.

Other Characteristics:

- Ideal for sterile water storage, polypropylene (PP) carboys can be autoclaved before refilling to stop recurring bacterial growth.
- Excellent chemical resistance
- Working temperature 0 to 135 °C
- Autoclavable

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **December 2015**