



Product Specification Sheet

VWR HIGH PERFORMANCE 50ml FREESTANDING CENTRIFUGE TUBES

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Screw caps	Packed	Type	Pk
525-0160	flat	bulk (tubes caps in separate bags)	Non sterile	500
525-0224	flat	50/bag, 10 bags/case	Sterile	500
525-0417	plug style	50/bag, 10 bags/case	Non sterile	500
525-0418	plug style	50/bag, 10 bags/case	Sterile	500
525-0416	plug style	bulk (tubes caps in separate bags)	Non sterile	500

Material: Polypropylene

Product description : molded graduations and white solvent resistant writing area

Product Dimensions: Length of Tube: **115 mm**
Outside Diameter: **29mm**
Inside Diameter: **25mm**

Pyrogens: certified testing performed lot by lot; Tests include Nuclease (RNase/DNase) and Endotoxin

Performance Testing : Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **April 2013**

VWR International bvba/sprl, Researchpark Haasrode 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium
<http://www.vwr.com>



Version: 2.0